


Evolution.

Revolution.



Innovated for you — to leave tape behind.

The AnchorFast SlimFit™ oral endotracheal tube fastener provides a convenient means to hold an oral endotracheal (ET) tube securely in place without the use of adhesive tape and is designed with the intent to accommodate smaller faces.

AnchorFast SlimFit
Oral Endotracheal Tube Fastener

AnchorFast SlimFit

Oral Endotracheal Tube Fastener

NEW
designed to
accommodate
smaller faces

NEW
upper lip
foam

Upper lip foam
designed to increase
surface area contact

Gliding tube shuttle
allows side-to-side
repositioning of the
ET tube for oral care,
while the tube remains
secure

Thumb pad
with increased
surface area for
finger or thumb
placement
when opening

NEW
thumb
pad

Gray track
that differentiates
from other
AnchorFast
products

NEW
gray
color

Skin barrier pads
primary means to secure the
device to the patient's face

Adjustable neck strap device
secondary securement designed
for adjustment at an accessible
area on the patient's head

NEW
neck
strap

Product/Description	Stock No	Quantity/Box
AnchorFast SlimFit oral endotracheal tube fastener Secures ET Tubes ranging from 5.0–8.0 mm	9787	12

For more information or to request a sample, please visit www.hollister.co.uk/anchorfast

For training requests, please contact Shekhar Nevgi on
Mob: 07881 787 885, Email: shekhar.nevgi@bywayconsulting.com



Hollister Limited
Rectory Court
42 Broad Street
Wokingham
Berkshire RG40 1AB

For more information and to
request FREE samples:
Tel: 0800 521 377
Email: samples.uk@hollister.com

For placing orders:
Tel: 01189 895088
Fax: 01189 895089
Email: hollister.tradeuk@hollister.com
www.hollister.co.uk

