



Discretion by Design

Introducing a discreet catheter with design that is truly inspired.



Infyna Chic

Intermittent Catheter



Continence Care

The art of discretion

A unique combination of colour, design, and discretion to help blend in with her life.

The Infyna Chic catheter was designed with input from clinicians to help provide a high level of discretion for women who use catheters. It is beautiful and easy to use. The Infyna Chic catheter can help a woman feel better about having to use one.



actual size




Easy to open with just one hand



Easy to grip and remove from case



Does not leak when reclosed



Designed with her in mind

Features that help change her intermittent catheter experience.



- Flip-top cap can be opened and closed with one hand
- Ready to use, hydrophilic catheter
- Just right stiffness helps support easy, touch-free insertion
- Catheter length helps give confidence that the bladder is fully drained
- Can attach a collection bag without removing the catheter from the case

Modern, feminine, discreet





Beautiful packaging and funnels in pastel colours—it's the details that make the difference.





Infyna Chic

Intermittent Catheter

Size (Charriere)	Catheter length	Tip	SKU	Quantity per box
 8 Ch	13 cm	Nelaton	7008	30
 10 Ch	13 cm	Nelaton	7010	30
 12 Ch	13 cm	Nelaton	7012	30
 14 Ch	13 cm	Nelaton	7014	30



actual size



To learn more or to request a sample, visit Hollister.co.uk/InfynaChic

Prior to use, be sure to read the Instructions For Use for information regarding Intended Use, Contraindications, Warnings, Precautions, and Instructions.

The Hollister logo and Infyna Chic are trademarks of Hollister Incorporated. ©2018 Hollister Incorporated. CC0551

